



**PRODUCT CATALOG**  
MEDICAL  
COMPRESSION STOCKINGS





## ONLINE SHOP

EXPERTS.BAUERFEIND.COM



## SERVICE GUARANTEES

### INTERNATIONAL CUSTOMER SERVICE



#### 100% AVAILABILITY

By telephone within our service times  
(Mo. – Th. 7.00 am – 05.00 pm, Fr. 7.00 am – 04.00 pm CET).  
By fax and email 24 hours per day.



#### ELECTRONIC ORDERING SYSTEM

Make use of our offers for a fast, faultless and effort  
reduced data transfer: online shop, wheel IT.



#### DESPATCH HANDLING

At your request we take care of the full processing  
of your shipment (incl. shipment quotations, additional  
documents) and we can bundle your orders flexible.



#### PROCESSING OF RETURNS

Within 5 days after receipt of goods.



#### ONE FACE TO THE CUSTOMER

You have one fixed contact person for all of  
your desires and concerns.



#### PROACTIVE INFORMATION

We keep you posted – you receive from us  
all important information about your order  
and delivery status.



#### RETURN OF GOODS

Redemption of faultless goods within  
30 days after delivery with 100 % refund  
or replacement (prior notification,  
filled form of return is required).

# CONTENTS

<b>4</b>	<b>About Bauerfeind</b>
<b>6</b>	<b>The four pillars of compression therapy</b>
<b>7</b>	<b>Systematic approach to compression</b>
<b>8</b>	<b>The new color system</b>
<b>10</b>	<b>Edging and top bands</b>
<b>12</b>	<b>Lymphology</b>
<b>14</b>	<b>Bauerfeind measurement technology</b>
<b>16</b>	<b>PREVENTION</b>
18	<b>VenoTrain® act</b> sheer elegance
20	<b>VenoTrain® act</b> cotton
20	<b>VenoTrain® act</b> microfiber
<b>22</b>	<b>STANDARD PHLEBOLOGICAL TREATMENT</b>
23	Perfect-Fit size chart
24	<b>VenoTrain® discrétion</b>
26	<b>VenoTrain® look</b>
28	<b>VenoTrain® micro</b>
30	<b>VenoTrain® micro</b> Design Edition and Batik
32	<b>VenoTrain® business</b>
34	<b>VenoTrain® soft</b>
36	<b>VenoTrain® pure</b>
36	VenoTrain pure size chart
<b>38</b>	<b>SPECIAL PHLEBOLOGICAL TREATMENT</b>
40	<b>VenoTrain® soft S</b>
40	Perfect-Fit size chart
42	<b>VenoTrain® cocoon</b>
42	Perfect-Fit size chart
46	<b>VenoTrain® impuls</b>
46	V6+ size chart
48	<b>VenoTrain® ulcertec</b>
48	X5+ size chart
<b>50</b>	<b>LYMPHOLOGICAL TREATMENT</b>
	Overview*
<b>52</b>	<b>ACCESSORIES</b>
<b>56</b>	<b>INDIVIDUAL CUSTOMER SERVICE AND ORDER</b>
56	Circular knit product range
58	Additional services   Alterations   Repairs
60	Order form   Order pads
60	Terms and Conditions of Sale and Delivery

\* Detailed information about products for lymphedema and lipedema can be found in the separate Lymphology Product Catalog (item number 271 4210 1001 001) as well as in the current price list.

# BAUERFEIND. BY PEOPLE FOR PEOPLE. FOR A BETTER QUALITY OF LIFE.

Bauerfeind AG is one of the world's leading manufacturers of medical aids. With our products and treatment concepts, we enable people to be more active and thus make an important contribution toward maintaining and restoring health.

## MEDICAL KNOW-HOW AND TECHNICAL EXPERTISE

Bauerfeind combines decades of expertise in technology and manufacturing with innovation: the tradition of a medium-sized, family-run business with a global mindset. Bauerfeind's specialists have worked very closely with physicians, scientists, and developers right from the start. The result: patented, award-winning products and services.

From the thread to the completed product, our compression stockings are developed and manufactured in Germany.


## EFFECTIVE VARIETY

Our VenoTrain compression stockings cover all uses/indications in modern compression therapy, offering medically proven effectiveness and unique material properties.


Thanks to the wide range of colors available, the compression products blend seamlessly into their users' everyday lives and can be combined with various fashion styles. For prevention or the treatment of lymphedema: the indication-specific solutions include a wide variety of circular knit and flat knit designs that are available as custom-made or standard products.

## QUALITY WITHOUT COMPROMISE

Bauerfeind products don't simply meet standards – they fulfill the most stringent requirements in terms of quality, effectiveness, and wearing comfort. This is confirmed by regular inspection procedures, for which our products are permitted to bear the seals and marks after successful testing. All VenoTrain products conform to the requirements of the German Medical Devices Act (MPG) and bear the RAL quality mark. As a seal of quality, the STANDARD 100 by OEKO-TEX® guarantees testing for harmful substances.







**OEKO-TEX®**  
 CONFIDENCE IN TEXTILES  
**STANDARD 100**  
 598-0153 HOHENSTEIN HTTI

Tested for harmful substances,  
[www.oeko-tex.com/standard100](http://www.oeko-tex.com/standard100)

**GÜTEZEICHEN**  
**MEDIZINISCHE**  
**KOMPRESSIONS-**  
**STRÜMPFE**  
 RAL





# THE FOUR PILLARS OF COMPRESSION THERAPY

## TREATMENT OPTIONS BASED ON MEDICAL INDICATIONS



### 01 PREVENTION (page 16–21)

The VenoTrain act **support stockings** offer beneficial relief for the legs and prevent vein problems. They improve circulation and keep legs fit when sitting or standing for long periods of time.



### 02 STANDARD PHLEBOLOGICAL TREATMENT (page 22–37)

The standard phlebological treatment is suitable for **mild to moderate venous disorders**. It also includes situations in which there is a particular risk for the development of vein problems, such as during pregnancy or traveling.

VenoTrain compression stockings, which have been developed precisely for these cases, improve venous activity and combat venous leg problems.



### 03 SPECIAL PHLEBOLOGICAL TREATMENT (page 38–49)

Some patients with vein problems have **special requirements** that call for compression stockings with unique modes of action and properties. These requirements can be very diverse.

That is why the VenoTrain special case range offers compression stockings especially for patients with restricted mobility, dry skin, pronounced venous insufficiency or venous leg ulcers.




### 04 LYMPHOLOGICAL TREATMENT (page 50–51)

Lymphedema or lipedema patients are generally dependent on compression therapy for their entire lives. The effectiveness and wearing comfort of the stockings are therefore particularly important, significantly impacting the patient's quality of life. VenoTrain flat knit and circular knit products generate the required therapeutic compression, whilst offering a unique wearing comfort, and help the patient live an active life.

# A SYSTEMATIC APPROACH TO COMPRESSION FOR THE RIGHT CHOICE

The Bauerfeind systematic approach will help you find the right VenoTrain compression stocking for your customers. The right choice depends on the severity of the condition.



\* Partsch, H.: Klassifikation von chronischen Venenerkrankungen heute, Vasomed 1 / 2003

\*\* Rabe, E. et al.: Bonner Venenstudie der DGP, Phlebologie 1 / 2003







# THE NEW COLOR SYSTEM

## SIMPLE. CLEAR. COORDINATED.

### SELECTING THE RIGHT COLOR HAS NEVER BEEN EASIER

The new VenoTrain color system makes it even easier for you and your customers to select the right color: Now, the circular knit VenoTrain medical compression stockings for mild to advanced vein problems – VenoTrain look, VenoTrain micro, VenoTrain soft, VenoTrain business and VenoTrain impuls – as well as the flat knit product VenoTrain curaflow, are all available in the colors Caramel, Navy, Anthracite and Black.

The popular variants VenoTrain micro and VenoTrain soft are also offered in the light skin tone Cream, the darker skin tone Amaretto, and the warm brown tone Espresso.

VenoTrain®	STANDARD COLORS				ADDITIONAL COLORS		
	Caramel	Navy	Anthracite	Black	Cream	Amaretto	Espresso
							
look	●	●	●	●			
micro	●	●	●	●	●	●	●
business	●	●	●	●			
soft	●	●	●	●	●	●	●
impuls	●	●	●	●			
curaflow	●	●	●	●	●		

### THE BENEFITS OF THE NEW COLOR SYSTEM FOR YOU:

- Coordinated color selection for the entire VenoTrain range
- Easy and quick overview of all products and their corresponding colors
- Personalized and straightforward customer advice

# EDGING AND TOP BANDS

## PERFECT HOLD WITH A HIGH LEVEL OF COMFORT

Our edging and top bands make sure that compression stockings fit securely and don't slip. Their soft materials are breathable, gentle on the skin and very comfortable to wear. Learn more about the individual styles and different designs.

### EDGING

#### SECURE AND COMFORTABLE



#### TOP BAND

For all knee-high compression stockings

The soft edging ensures wearing comfort and doesn't dig into the skin. Its neutral look works for both women and men.



#### COMFORT BANDS

For the Venotrain look and Venotrain soft pantyhose (standard)

This edging provides an excellent fit and a high level of comfort thanks to its higher elasticity and wide band. It sits comfortably without constricting.

### SILICONE TOP BANDS

#### SECURE AND INDIVIDUALIZED



#### SILICONE MICRO DOTS

For Venotrain look (standard and custom-made)

The highly elastic silicone top band with its small, flat silicone dots is very inconspicuous, even when worn under tight clothing, and leaves no marks on the skin. The skin-friendly material with its soft surface adapts perfectly to the leg contour and is very comfortable to wear even with soft connective tissue. The silicone top band comes in one width for each size and no longer needs to be produced in both the regular and wide size.



#### SILICONE DOTS

For Venotrain micro, Venotrain soft, Venotrain pure, Venotrain soft S and Venotrain impuls (standard and custom-made)

Our classic silicone top band is mainly used for thigh-high stockings. It is gentle on the skin and ideal for men thanks to its neutral design.





### GARTER BELT

For Venotrain soft and Venotrain impuls (standard and custom-made)

The garter belt without silicone coating improves the hold of the stockings and offers plenty of freedom of movement. It is suitable for sensitive skin and can be worn on one or both sides.



### SILICONE DOT PATTERN

For Venotrain look (custom-made), Venotrain micro, Venotrain soft and Venotrain soft S (standard and custom-made)

This skin-friendly silicone top band with its elegant pattern is especially popular with women.



### LACE

For Venotrain soft (custom-made)

The curved, flat coating is breathable and provides a good hold. Its lace design is especially appealing to women.



### LACE SENSITIVE

For Venotrain look and Venotrain micro (standard and custom-made)

The highly elastic silicone top band has an ultra-thin, continuous coating that provides a perfect hold. It adapts perfectly to the leg contour, even with soft connective tissue. The silicone top band comes in one width for each size and no longer needs to be produced in both the regular and wide size.

# HELLO NEW ATTITUDE TOWARD LIFE!

## A NEW LOOK FOR LYMPHOLOGICAL THERAPY

Every day is more pleasant when you're feeling better. That is why Bauerfeind tirelessly works on product innovation.



### + NEW TOE CAP, NOW ALSO WITHOUT IRRITATING SEAMS

- + GLOVE WITHOUT IRRITATING FINGER SEAMS
- + ARM SLEEVE WITH AN EXTRA-SOFT INNER ELBOW AREA AND ANATOMICAL ELBOW
- + SPECIALLY DESIGNED COMFORT ZONES
- + EXTRA SOFT BEHIND THE KNEE
- + SEAMLESS AND ANATOMICALLY SHAPED 3D CHEST AREA, TOE WITH SOFT TIP AS STANDARD


### FOR MORE MOBILITY AND COMFORT DURING EVERYDAY ACTIVITIES

With the VenoTrain curaflow, we are offering flat knit products that ensure unique wearing comfort when therapeutic compression is required. For patients with lipedema and lymphedema, our products make everyday activities less complicated.


### STYLISH CONTRASTS

The different colors are very popular with women who wear compression stockings because they can be combined with a variety of outfits. Anthracite, in particular, represents classic elegance. Exciting contrasts can be created with the bordeaux hue.

Detailed information about products for lymphedema and lipedema can be found in the Lymphology Product Catalog (item number 271 4210 1001 001) as well as in the current price list.







# QUALITY THAT CAN BE MEASURED.

## BAUERFEIND MEASUREMENT TECHNOLOGY



### BODYTRONIC® 200

#### LEG HEALTH CAN BE MEASURED

- Early detection of vein problems that cannot be seen from the outside
- A quick, non-invasive vein function measurement process using infrared sensors



### BODYTRONIC® 410

#### FLEXIBLE 3D SCANNER FOR THE LOWER EXTREMITIES

- Contact-free measurement for compression stockings, with no need for special measurement clothing
- 3D image of the measurement area up to the waist
- Quick and easy to set up and take down



### BODYTRONIC® 400

#### STATIONARY 3D SCANNER FOR THE LOWER EXTREMITIES INCLUDING FOOT SCANNER

- Contact-free measurement for compression stockings, with no need for special measurement clothing
- 3D image of the measurement area up to the waist
- System available with or without foot scanner (for foot orthoses)






 NEW


**BODYTRONIC® 610**  
 FULL-BODY SCANNER  
 INCLUDING PRODUCT SUPPLY  
 FOR ALL EXTREMITIES

- Contact-free measurements for medical compression stockings, knee and foot supports, knee orthoses as well as for back, arm and shoulder products without the need for special measurement clothing
- Optional full-body scan, export of 3D data
- Integrated leg volume measurement
- Standing space of only 1.60 m<sup>2</sup>

**01**  
 FUNCTIONAL  
 MEASUREMENT  
 SYSTEMS


**02**  
 PRECISE  
 MEASURE-  
 MENT

**03**  
 PERFECTLY  
 FITTING  
 PRODUCTS



**BODYTRONIC® 600**  
 FULL-BODY SCANNER INCLUDING  
 PRODUCT SUPPLY FOR THE LOWER  
 EXTREMITIES

- Contact-free measurement for compression stockings, knee and foot supports as well as knee orthoses, with no need for special measurement clothing
- Optional upper body scan, export of 3D data
- Integrated leg volume measurement
- Standing space of 2.80 m<sup>2</sup>



# PREVENTION

## GENTLE SUPPORT AT WORK AND LEISURE

Support stockings are an effective means against tired and heavy legs. They provide excellent strain relief and counteract impaired circulation. New skin-friendly materials meet stylish designs and a wide choice of colors.













# VenoTrain<sup>®</sup> act sheer elegance








## PROVIDING ELEGANT SUPPORT

In 70 or 140 denier, the VenoTrain act sheer elegance provides light to medium support in a modern translucent design that is available in a wide range of colors and styles.

### FEATURES

- Elegant and transparent
- Available in seven shades
- Available as a knee-high or thigh-high stocking with a silicone lace band, and as pantyhose
- Also available as pantyhose for pregnant women with a special panty area that expands as the mother's bump grows
- Calibrato (CAL) pantyhose for particularly feminine curves

### STANDARD

		AD	AG	AT	AT CAL	ATU
<b>SIZES</b>		1-4	1-5	1-5	2-5	1-5
<b>DENIER</b>		70/140	70/140	70/140	70/140	140
<b>COLORS</b>	 SAND	●	●	●		
	 CAMEL	●	●	●	●	●
	 AMARETTO	●	●	●		
	 ESPRESSO	●	●	●		
	 BLACK	●	●	●	●	●
	 GRAPHITE	●	●	●		
	 NAVY	●	●	●		

### DETAILS

TOE closed

#### MATERIAL CONTENT\*

**AG (70 DEN)** 81 % polyamide  
19 % elastane


**AG (140 DEN)** 73 % polyamide  
27 % elastane

#### CARE INSTRUCTIONS




### SIZE CHARTS

#### AG, AT, ATU



#### AT CAL




#### AD

#### SHOE SIZES

1	37-39
2	39-41
3	41-43
4	43-45

\* Additional material contents are available on request.

No alterations, repairs, or additional services are offered for the VenoTrain act. It cannot be returned in exchange for a credit note or another item.





COTTON



MICROFIBER



COTTON



COTTON



# VenoTrain<sup>®</sup> act cotton, microfiber

## PROVIDING GENTLE SUPPORT

VenoTrain act cotton and VenoTrain act microfiber, each available for women and men, offer a high degree of support in 280 denier and are reliable companions when traveling, at work or leisure.

### DETAILS

**TOE** closed

#### MATERIAL CONTENT

- COTTON: 27 % polyamide, 10 % elastane, 63 % cotton
- MICROFIBER: 94 % polyamide, 6 % elastane

#### CARE INSTRUCTIONS



### FEATURES

#### COTTON

- Particularly high level of wearing comfort based on cotton
- With a distinctive ribbed look
- Opaque
- Highly supportive knee-high stocking (280 denier)
- Available in four colors
- For men and women



#### MICROFIBER

- Particularly soft and gentle on the skin, thanks to the high microfiber content
- Fits snugly around the leg, while allowing the skin to breathe
- Highly supportive knee-high stocking (280 denier)
- Available in four colors
- For men and women

### STANDARD

		
		AD
<b>SIZES</b>		1 – 4
<b>DENIER</b>		280
<b>COLORS</b>	 NATURE	●
	 BLACK	●
	 GRAPHITE	●
	 NAVY	●

### SIZE CHART

		
	AD	AD
<b>SHOE SIZES</b>	<b>WOMAN</b>	<b>MAN</b>
1	37 – 39	39 – 41
2	39 – 41	41 – 43
3	41 – 43	43 – 45
4	43 – 45	45 – 47

No alterations, repairs, or additional services are offered for the VenoTrain act. It cannot be returned in exchange for a credit note or another item.





# TREATMENT IN PHLEBOLOGICAL STANDARD CASES

## FOR MILD TO MODERATE VENOUS PROBLEMS

The standard phlebological treatment is suitable **for mild to moderate venous disorders**. It also includes situations in which there is a particular risk for the development of vein problems, such as during pregnancy or traveling. VenoTrain compression stockings, which have been developed precisely for these cases, improve venous activity and combat venous leg problems.

## PERFECT-FIT

### STANDARD SIZE CHART FOR CIRCULAR KNIT PRODUCTS

VenoTrain® look, VenoTrain® micro, VenoTrain® business, VenoTrain® soft



	CIRCUMFERENCE	S		M		L		XL	
		REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS
<b>Length I/G</b> short 61–70 long 70–80 Type AG, AT	cG silicone top band, wide	51–57	59–66	56–63	62–69	60–67	66–74	63–71	71–80
	cG silicone top band, regular	45–51	52–59	49–56	55–62	53–60	58–66	55–63	62–71
	cG top band, AT, silicone top band*	45–57	52–66	49–63	55–69	53–67	58–74	55–71	62–80
<b>Length I/D</b> short 35–40 long 40–46 Type AD	cF	40–52	43–55	43–55	46–60	45–59	50–64	49–63	52–68
	cE	33–40	36–44	35–43	39–49	38–48	42–53	42–53	47–58
	cD	29–35	32–38	32–38	35–42	34–42	37–45	37–45	41–49
	cC	31–37	34–40	34–40	37–44	37–45	39–47	40–48	43–51
	cB1	25–29	27–31	27–32	29–34	29–35	31–37	32–38	34–40
	cB	20–22		22–25		25–28		28–31	
	cY	29–33		31–35		34–40		37–44	
	cA	19–24		22–27		25–31		28–35	

\* Lace "sensitive" silicone top band for VenoTrain look and VenoTrain micro  
Micro silicone dots top band for VenoTrain look



### STANDARD SIZE CHART: VenoTrain discrétion

#### LEG CIRCUMFERENCE

Size	1	2	3	4
<b>cG</b>	42-56	49-63	56-70	63-77
<b>cC</b>	28-36	32-40	35-43	38-46
<b>cB</b>	18-21	21-24	24-27	27-31

#### LENGTH

Size	short	long
<b>AG   AT</b>	< 70	> 70
<b>AT</b>	< 38	> 38



# VenoTrain® discrétion

## THE TRANSPARENT STOCKING

The exceptionally fine material of these compression stockings gives them a stylish transparency and lightness.

Available in various designs and versatile colors, these stockings will give you an extra sense of vitality and well-being in your work, leisure and everyday life.

### COLORS



### FEATURES

- Highly elastic, making the stockings easy to put on
- Shapes the legs without constricting them
- Breathable fine knit for a long-lasting wearing comfort

### DETAILS

**SIZING SYSTEM** VenoTrain discrétion (I–IV)

#### COMPRESSION CLASSES

- Class 1 (18–21 mmHg)
- Class 2 (23–32 mmHg)

**TOE** closed, open

#### MATERIAL CONTENT

- Class 1: 75 % polyamide, 25 % elastane
- Class 2: 70 % polyamide, 30 % elastane

#### CARE INSTRUCTIONS



### ● STANDARD

### ○ CUSTOM-MADE SOLUTIONS

	 <b>AD</b> Knee stockings	 <b>AG</b> Thigh-high stockings with silicone lace band	 <b>AG</b> Thigh-high stockings with silicone dotted band	 <b>AT</b> Pantyhose	
<b>EDGING</b>					
<b>SILICONE TOP BAND</b>					
Clair	●*	●*	●*		
Beige	●○	●○	●○	●○	
Espresso	●○	●○	●○	●○	
Black	●○	●○	●○	●○	
Java	●	●	●	●	

\* only available in classe 2 (15–20 mmHg)  
available in selected countries







# VenoTrain® look

## ATTRACTIVE AND GLOSSY

The VenoTrain look is a real highlight among medical compression stockings thanks to its glossy look – for a stylish appearance on any occasion!

The VenoTrain look enhances the figure with a push-up effect. Stretch comfort and its special trilobal fiber guarantee a perfect fit.

### STANDARD COLORS



### FEATURES

- Fashionable and elegant gloss effect
- Slightly transparent

### DETAILS



#### SIZING SYSTEM PERFECT-FIT (S–XL)

#### COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

#### TOE closed

#### MATERIAL CONTENT




- Ccl 1: 60 % polyamide, 40 % elastane
- Ccl 2: 55 % polyamide, 45 % elastane

#### CARE INSTRUCTIONS



#### ● STANDARD

#### ○ CUSTOM-MADE SOLUTIONS

EDGING	 AD	 AG	 AT
<b>SILICONE TOP BAND</b>			
Top band	● ○		● ○
Lace "sensitive"		● ○	
Micro silicone dots		● ○	
Silicone dot pattern		○	





# VenoTrain® micro

## FASHIONABLE, MODERN AND VERSATILE

With a particularly high microfiber content that makes it comfortable to wear, the VenoTrain micro is one of the most popular compression stockings.

Available in many designs and a wide selection of colors, it highlights the individual style of your customers.

### STANDARD COLORS



### ADDITIONAL COLORS



### FEATURES

- More than 50 per cent microfiber content: breathable, skin-friendly, and easy-care
- Figure-shaping panty area (pantyhose) and comfort sole
- Opaque full-fashioned stocking quality

### DETAILS



#### SIZING SYSTEM PERFECT-FIT (S-XL)

#### COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

TOE closed, open

#### MATERIAL CONTENT

- Ccl 1: 65 % polyamide, 35 % elastane
- Ccl 2: 55 % polyamide, 45 % elastane

#### CARE INSTRUCTIONS



#### ● STANDARD

#### ○ CUSTOM-MADE SOLUTIONS

EDGING	AD	AG	AT	ATU
SILICONE TOP BAND				
Top band	● ○		● ○	● ○
Lace "sensitive"		● ○		
Silicone dots	● ○	● ○		
Silicone dot pattern		● ○		



Select favorite designs quickly and easily:

[www.bauerfeind.com/venotrain-stylist](http://www.bauerfeind.com/venotrain-stylist)



Select favorite designs  
quickly and easily:  
[www.bauerfeind.com/venotrain-stylist](http://www.bauerfeind.com/venotrain-stylist)







TANGO

# VenoTrain® micro DESIGN EDITION AND BATIK

## THE SECOND SKIN FOR TRENDSETTERS

A true fashion statement: the Design Edition and individually dyed VenoTrain micro Batik versions create real highlights.


With a subtle Tango or Jive pattern, or wildly textured in batik with navy or black, your customers' legs will be real eye-catchers.



BATIK NAVY



JIVE



BATIK BLACK

### STANDARD COLORS

	TANGO*	JIVE*	BATIK NAVY	BATIK BLACK
CARAMEL				
ANTHRACITE				
NAVY				
BLACK				

### ADDITIONAL COLORS

	TANGO*	JIVE*	BATIK NAVY	BATIK BLACK
CREAM				
AMARETTO				
ESPRESSO				

- STANDARD
- CUSTOM-MADE SOLUTIONS

<b>EDGING</b>				
<b>SILICONE TOP BAND</b>	<b>AD</b>	<b>AG</b>	<b>AT</b>	<b>ATU**</b>
Top band	● ○		● ○	● ○
Lace "sensitive"		● ○		
Silicone dots**	● ○	● ○		
Silicone dot pattern		● ○		

\* The Design Edition is available only with closed toes.

\*\* Not available as Design Edition





# VenoTrain® business

## TIMELESSLY ELEGANT

While traveling or at work: the VenoTrain business meets every requirement for a premium knee-high stocking. The combination of cotton and Timbrelle® high-tech fiber offers top-class wearing comfort with thermal regulation. Highly elastic LYCRA® and the extra-wide supportive band ensure that the stocking is easy to put on and stays perfectly in place.

### FEATURES

- Innovative mix of materials
- Reinforced sole
- Timbrelle® high-tech fiber providing long-lasting freshness and healthy foot conditions



### STANDARD COLORS



### DETAILS



#### SIZING SYSTEM PERFECT-FIT (S–XL)

#### COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

#### TOE closed

#### MATERIAL CONTENT

- Ccl 1: 80 % polyamide, 14 % elastane, 6 % cotton
- Ccl 2: 71 % polyamide, 24 % elastane, 5 % cotton

#### CARE INSTRUCTIONS



- **STANDARD**
- **CUSTOM-MADE SOLUTIONS**

#### EDGING SILICONE TOP BAND



AD

Top band



Silicone dots



Available in shoe sizes 36–41 (foot short) or 41–46 (foot long)





# VenoTrain® soft

## GENTLE AND EFFECTIVE

The VenoTrain soft, always the right choice: this versatile stocking can withstand plenty of wear, yet it is very soft and gentle on the skin.

It is available in many designs with a proven medical effect and a beneficial massage action.

### STANDARD COLORS



### ADDITIONAL COLORS



### FEATURES

- For men and women
- Particularly gentle on the skin
- Also available as a single stocking with a standard design in caramel
- Also available as pantyhose for pregnant women with a special panty area that expands as the mother's bump grows (ATU)

### DETAILS



**SIZING SYSTEM** PERFECT-FIT (S–XL)

#### COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

**TOE** closed, open, no foot (custom-made), stump (custom-made)

#### MATERIAL CONTENT

- Ccl 1: 70 % polyamide, 30 % elastane
- Ccl 2: 65 % polyamide, 35 % elastane

#### CARE INSTRUCTIONS



### ● STANDARD

### ○ CUSTOM-MADE SOLUTIONS

EDGING	AD	AG	AGHB	AT	ATU	ATH	ATX	ATE
SILICONE TOP BAND	● ○			● ○	○	○	○	○
Top band								
Lace		○						
Silicone dots	● ○	● ○						
Silicone dot pattern		● ○						
Garter belt			● ○					

C0

C1

C2

C3

C4

C5

C6



The VenoTrain pure is available as a knee-high or thigh-high stocking with an open or closed toe.

### STANDARD SIZE CHART VenoTrain® pure

	CIRCUMFERENCE	S		M		L		XL	
		REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS
<b>Length lG</b>	<b>cG</b>	46-55	55-64	49-58	58-67	52-62	62-72	55-66	66-78
<b>short</b> 61-70									
<b>long</b> 70-80									
<b>Type</b> AG	<b>cE</b>	33-40	36-44	35-43	39-49	38-48	42-53	42-53	47-58
<b>Length lD</b>	<b>cD</b>	29-38		32-42		34-45		37-49	
<b>short</b> 35-40	<b>cC</b>	31-40		34-44		37-47		40-51	
<b>long</b> 40-46									
<b>Type</b> AD	<b>cB</b>	20-22		22-25		25-28		28-31	
	<b>cY</b>	29-33		31-35		34-40		37-44	
	<b>cA</b>	19-24		22-27		25-31		28-35	



# VenoTrain® pure

## BASIC CARE

For uncompromising medical effectiveness with a simplified sizing system providing quick and economical treatment. With a low-key natural color, it is suitable for him and her as well as gentle on the skin and comfortable to wear.

### COLOR

NATURAL



### FEATURES

- Soft, durable material
- Available in four standard sizes, each in two lengths
- Also available as a single stocking

### DETAILS

**SIZING SYSTEM** VenoTrain pure (S–XL)

#### COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

**TOE** closed, open

#### MATERIAL CONTENT

- Ccl 1: 77 % polyamide, 23 % elastane
- Ccl 2: 74 % polyamide, 26 % elastane

#### CARE INSTRUCTIONS



### ● STANDARD

#### EDGING SILICONE TOP BAND



AD



AG

Top band



Silicone dots



The VenoTrain pure cannot be returned in exchange for a credit note or another item.

C0

C1

C2

C3

C4

C5

C6





# TREATMENT IN PHLEBOLOGICAL SPECIAL CASES FOR MODERATE TO SEVERE VEIN PROBLEMS

Some patients with vein problems have **special requirements** that call for compression stockings with unique modes of action and properties. These requirements can be very diverse. That is why the Venotrain special range offers compression stockings especially for patients with:

- **RESTRICTED MOBILITY**  
**Venotrain soft S**
- **DRY SKIN**  
**Venotrain cocoon**
- **PRONOUNCED VENOUS INSUFFICIENCY**  
**Venotrain impuls**
- **VENOUS LEG ULCERS**  
**Venotrain ulcertec**



The VenoTrain soft S is available as a knee-high or thigh-high stocking with an open or closed toe.

### PERFECT-FIT STANDARD SIZE CHART

	CIRCUMFERENCE	XS*		S		M		L		XL	
		REGULAR	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	
<b>Length lG</b>											
<b>short</b>	61-70										
<b>long</b>	70-80										
<b>Type</b>	AG, AT										
		<b>cG</b> silicone top band, wide		51-57	59-66	56-63	62-69	60-67	66-74	63-71	71-80
		<b>cG</b> silicone top band, regular	44-50	45-51	52-59	49-56	55-62	53-60	58-66	55-63	62-71
		<b>cF</b>	39-50	40-52	43-55	43-55	46-60	45-59	50-64	49-63	52-68
		<b>cE</b>	32-38	33-40	36-44	35-43	39-49	38-48	42-53	42-53	47-58
		<b>cD</b>	28-31	29-35	32-38	32-38	35-42	34-42	37-45	37-45	41-49
		<b>cC</b>	29-33	31-37	34-40	34-40	37-44	37-45	39-47	40-48	43-51
		<b>cB1</b>	23-26	25-29	27-31	27-32	29-34	29-35	31-37	32-38	34-40
		<b>cB</b>	18-20	20-22		22-25		25-28		28-31	
		<b>cY</b>	27-31	29-33		31-35		34-40		37-44	
		<b>cA</b>	17-22	19-24		22-27		25-31		28-35	

\* only for VenoTrain soft S – limited product range



# VenoTrain® soft S

## COMFORT IN THE EVENT OF RESTRICTED MOBILITY

Based on proven VenoTrain soft quality which offers a gentle massage effect to stimulate the blood flow, the VenoTrain soft S also takes into account special requirements of people who spend most of their time seated.

The wide heel and comfortable instep area ensure that the foot, which is angled during sitting, has enough space, and that no creases can develop.

The robust knitted fabric can withstand plenty of wear, and can also be easily put on by care staff.

### FEATURES

- Wide heel area
- Adapted to the typical anatomical foot posture while sitting (90-degree angle)
- Breathable material with a revitalizing effect

### STANDARD COLORS

Standard products only available in caramel

CARAMEL 

BLACK 

NAVY 

### ADDITIONAL COLORS

CREAM 

ESPRESSO 



### DETAILS

**SIZING SYSTEM** PERFECT-FIT (XS – XL)

#### COMPRESSION CLASSES

- Ccl 1 (18–21 mmHg)
- Ccl 2 (23–32 mmHg)

**TOE** closed, open

#### MATERIAL CONTENT

- Ccl 1: 70 % polyamide, 30 % elastane
- Ccl 2: 65 % polyamide, 35 % elastane

#### CARE INSTRUCTIONS



- **STANDARD**
- **CUSTOM-MADE SOLUTIONS**

### EDGING

Top band

Silicone dots

Silicone dot pattern



AD



AG

● ○

● ○

● ○

● ○




SKIN REPAIR EFFECT



CLINICALLY PROVEN SKIN-FRIENDLINESS

PERFECT-FIT STANDARD SIZE CHART

	CIRCUMFERENCE	S		M		L		XL	
		REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS
<b>Length (D)</b>									
<b>short</b>	35-40								
<b>long</b>	40-46								
<b>Type</b>	AD								
<b>cD</b>		29-35	32-38	32-38	35-42	34-42	37-45	37-45	41-49
<b>cC</b>		31-37	34-40	34-40	37-44	37-45	39-47	40-48	43-51
<b>cB1</b>		25-29	27-31	27-32	29-34	29-35	31-37	32-38	34-40
<b>cB</b>		20-22		22-25		25-28		28-31	
<b>cY</b>		29-33		31-35		34-40		37-44	
<b>cA</b>		19-24		22-27		25-31		28-35	



# VenoTrain® cocoon

## FOR DRY AND SENSITIVE SKIN

The medical compression stocking with innovative dual function protects and soothes dry and sensitive skin.

High-quality skincare substances repair the natural skin barrier, making dry skin more supple.

### FEATURES

- Clinically tested skincare substances with cosmetic effect
- Clinically proven skin-friendliness
- With an innovative mix of natural fibers
- The lipid care complex improves skin hydration and counteracts itching and scaling
- For men and women

### COLORS



### DETAILS

**SIZING SYSTEM** PERFECT-FIT (S–XL)

**COMPRESSION CLASS**

- Ccl 2 (23–32 mmHg)

**TOE** closed

**MATERIAL CONTENT**

- Ccl 2: 52 % polyamide, 26 % elastane, 11 % cotton, 11 % cellulose

**CARE INSTRUCTIONS**



### ● STANDARD

**EDGING**  
**SILICONE TOP BAND**



**AD**

Top band



Available in shoe sizes 36–41 (foot short) or 41–46 (foot long)



# VenoTrain® cocoon

## THE COMPRESSION STOCKING WITH A DUAL EFFECT


It is common for patients wearing compression stockings to be affected by dry skin. Bauerfeind was able to develop a compression stocking which combines compression with long-lasting skincare, turning a compression stocking into a feel-good compression stocking. This is possible with the "magical" fibers contained in the VenoTrain cocoon.

**NEW**  
in two foot lengths


The VenoTrain cocoon medical compression stocking has a unique dual function for people with vein problems: compression and long-lasting skincare in one.

- + **With integrated cocoon skincare fibers made from cellulose, cotton, and a lipid skincare complex**
- + **Aids skin regeneration, preventing itching and scaling**


The care stocking is perfume-free and tested in accordance with Standard 100 by OEKO-TEX® as a guarantee of its safety.



59%  
of patients suffer from dry skin\*



33%  
of patients suffer from itching\*



84%  
of patients think it is important that compression stockings do not cause skin problems\*\*




## "WITH COCOON, CARE STARTS IN THE FIBERS"

"The stocking is called cocoon because we wanted to connect the properties with a word. And I think we managed to do that. A 'surrounding and protecting cover' fits this stocking perfectly," Gerti Stegmann explains, Product Developer in the Bauerfeind Innovation Center.

\* Int Angiol. 2009 Oct; 28(5): 385–93. Quality of life and patients view of compression therapy. Reich-Schupke S., Murmann F., Altmeyer P., Stücker M., Source Department of Dermatology, Ruhr University Bochum, Bochum, Germany, PMID: 19935593.

\*\* Allensbacher Archiv, IFD-Umfrage 9691

## THE TWO-IN-ONE MODE OF ACTION: SUPPORTS LONG-LASTING, NATURAL SKIN REGENERATION




### The Venotrain cocoon consists of a new mix of natural fibers with an innovative cocoon care thread.

This thread contains special lipid depots with high-quality care substances, which are continuously applied to the skin when the stocking is worn.

The clinically tested lipid care complex works directly on the outermost layer of the skin – the stratum corneum. It supports long-lasting, natural skin regeneration and effectively prevents itching and scaling. Compression and genuine skincare: the combination that makes the Venotrain cocoon new and unique.



Dry, sensitive skin  
before wearing the  
Venotrain cocoon



Smooth, healthy skin  
after repeatedly wearing  
the Venotrain cocoon

## "THIS COMBINATION IS COMPLETELY NEW!"

"The combination of materials, together with the active substance, is a completely new technology, never seen before."

Gerti Stegmann, Product Developer in the  
Bauerfeind Innovation Center

SKIN REPAIR  
EFFECT







V6+ STANDARD SIZE CHART

	CIRCUM-FERENCE	I		II		III		IV		V		VI	
		REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS
<b>Length ℓG</b>	<b>cG</b>	45-53	51-59	48-56	54-62	51-59	57-66	54-63	60-69	57-66	63-73	60-70	66-76
<b>short</b> 62-70	<b>cF</b>	39-45	43-50	42-49	46-53	45-52	49-57	48-56	52-60	51-59	55-64	54-63	58-67
<b>long</b> 70-80	<b>cE</b>	32-37	35-40	34-39	37-42	36-41	39-44	38-44	41-47	40-46	43-49	42-48	45-51
<b>Type</b> AG	<b>cD</b>	28-32	30-34	30-34	32-37	32-37	34-39	34-39	36-41	36-41	38-43	38-43	40-46
<b>Length ℓD</b>	<b>cC</b>	29-33	31-35	31-35	33-38	33-38	35-40	35-40	37-42	37-42	39-44	39-44	41-47
<b>short</b> 35-40	<b>cB1</b>	23-27	24-28	25-29	26-30	27-31	28-32	29-34	30-35	31-36	32-37	33-38	34-39
<b>long</b> 40-44	<b>cB</b>	19-20,5	20,5-22	20,5-22	22-24	22-24	24-26	24-26	26-28,5	26-28,5	28,5-31	28,5-31	
<b>Type</b> AD	<b>cY</b>	25-29	26-31	26-31	28-33	28-33	30-35	30-35	32-38	32-38	35-41	35-41	
	<b>cA</b>	19-22	20-23	20-23	22-25	22-25	24-27	24-27	26-30	26-30	28-33	28-33	




# VenoTrain® impuls

## POWERFUL SUPPORT FOR PRONOUNCED VENOUS INSUFFICIENCY

Produces a particularly deep action by exerting high working pressure with even the slightest muscle movements. Easy to put on and breathable thanks to microfibers.

Highly effective for advanced vein problems, post-thrombotic syndrome, follow-up care for healed ulcers, and therapy in the early stages of lymphedema (0–1).

### STANDARD COLORS



### FEATURES

- The short-stretch knitted fabric generates a high working pressure with even the slightest muscle movements
- Skin-friendly microfiber material makes the product breathable and ensures climate control
- Highly differentiated V6+ sizing system (I–VI)

### DETAILS

**SIZING SYSTEM** V6+ (I–VI)

**COMPRESSION CLASSES**

- Ccl 2 (23–32 mmHg)
- Ccl 3 (34–46 mmHg)

**TOE** closed, open

**MATERIAL CONTENT**

- Ccl 2: 80 % polyamide, 20 % elastane
- Ccl 3: 75 % polyamide, 25 % elastane

**CARE INSTRUCTIONS**



● **STANDARD**

○ **CUSTOM-MADE SOLUTIONS**


	AD	AG	AGHB	AT	ATH	ATX
<b>EDGING</b>						
<b>SILICONE TOP BAND</b>						
Top band	● ○			○	○	○
Silicone dots	● ○	● ○				
Garter belt			● ○			

C0	C1	C2	C3	C4	C5	C6
----	----	----	----	----	----	----



### X5+ STANDARD SIZE CHART

Length (D)	CIRCUM-FERENCE	XS		S		M		L		XL		
		REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	REGULAR	PLUS	
short 35-40	cC	24-29	28-34	27-32	31-38	29-36	35-42	33-40	38-47	36-44	43-52	
long 40-45		cB		20-21,5		21,5-23		23-25		25-28		28-32
Type AD												



# VenoTrain® ulcertec

## FOR VENOUS LEG ULCERS

The combination of an overstocking with a patented rhomboid knitted fabric and a skin-friendly liner for holding the wound dressing in place makes this compression stocking easy to put on without any need for assistance.

The liners are also ideal for providing a continuous low level of compression during the night. These are available to order separately in a pack of three and are suitable for washing at high temperatures.

### FEATURES

- Coordinated compression of overstocking and liner
- Unique compression profile that accelerates the healing process while ensuring exceptional patient comfort
- Differentiated treatment based on two designs within one compression class

#### 39 – moderate

- Compression class 3 (34–46 mmHg)
- Liner with gentle compression providing a compression pressure in the lower range of compression class 3 in combination with the overstocking

#### 46 – strong

- Compression class 3 (34–46 mmHg)
- Liner with moderate compression providing a compression pressure in the upper range of compression class 3 in combination with the overstocking

### COLORS

NATURAL  
OVERSTOCKING



WHITE  
LINER



### DETAILS

**SIZING SYSTEM** X5+ (XS–XL)

#### COMPRESSION CLASSES

- 39–moderate: 34–46 mmHg  
Liner with gentle compression
- 46–strong: 34–46 mmHg  
Liner with moderate compression

#### TOE

- Open (overstocking)
- Closed (liner)

#### MATERIAL CONTENT

##### LINER

- 39–moderate: 87 % polyamide, 13 % elastane
- 46–strong: 69 % polyamide, 31 % elastane

##### OVERSTOCKING

- 39–moderate / 46–strong: 79 % polyamide, 21 % elastane

#### CARE INSTRUCTIONS

- Overstocking:
- Liner:

#### ● STANDARD

#### ○ CUSTOM-MADE SOLUTIONS

**EDGING  
SILICONE TOP BAND**



AD

Top band







# LYMPHOLOGICAL TREATMENT

Bauerfeind has a wide range of compression products for lymphological conditions which take into account the patients' different requirements.



The circular knit **VenoTrain soft** arm sleeve is generally used to treat mild arm lymphedema, and to prevent imminent arm lymphedema after removal of the lymph node in the armpit or during radiation treatment.



The circular knit **VenoTrain impuls** compression stocking combines a wide range of benefits: its breathable microfiber and firm compression thread offer unique wearing comfort with a high degree of effectiveness. Only very little movement is needed for the stocking to exert the required therapeutic pressure.



The flat-knit **VenoTrain delight** compression stocking is indicated for use in cases of advanced venous disorders and lymphedema. The breathable double-stretch knit is skin-friendly and easy to put on and take off. Additionally, the firm material offers a good shaping effect and is very durable.



Thanks to our custom-made **VenoTrain curaflow** flat knit products, we can offer lymphatic patients perfectly tailored solutions with more than 1,000 different options of combining products and additional features. The beneficial therapeutic features of the range include the knitted fabric's high microfiber content, the flat sensitive stitching, and special comfort zones such as the extra-soft area behind the knee.

	LEG & FOOT	ARM & HAND	THORAX
<b>VenoTrain® soft</b> Arm		<b>Ccl 1 – 2</b>	
<b>VenoTrain® impuls</b>	<b>Ccl 2 – 3</b>		
<b>VenoTrain® delight</b>	<b>Ccl 2 – 4</b>		
<b>VenoTrain® curaflow</b>	<b>Ccl 1 – 4</b>	<b>Ccl 1 – 3</b>	<b>Ccl 1</b>

Comprehensive information and product details of our lymphology portfolio can be found in the Lymphology Product Catalog.

You are welcome to order it, indicating the item number 271 4210 1001 001, or download it directly using the following link:

[www.bauerfeind.com/downloads](http://www.bauerfeind.com/downloads)



A close-up photograph of a woman with dark hair, smiling warmly as she applies a white, textured cream to her right knee. She is wearing a white, ribbed, short-sleeved top. The background is softly blurred, showing what appears to be a window with light coming through. The overall mood is calm and caring.

## **ACCESSORIES**

### THE PERFECT ADDITION FOR EVERYDAY USE

These handy extras from the VenoTrain accessories range are useful for ensuring that your patients feel comfortable using and wearing their VenoTrain compression stockings. They provide effective treatment support, make it easier to handle the products, and care for the compression stockings as well as your customers' skin.





01



02



03



04



05



06

### 01 Lingerie bag

To protect the compression stockings during every machine wash

### 02 Bauerfeind care foam\*

Contains panthenol to protect and soothe particularly stressed and sensitive skin; before putting on the stockings and in the evening after taking them off

### 03 Detergent\*

Cleans gently, ensuring careful product treatment and providing lasting protection

### 04 Rubber gloves

Help with putting the compression stockings on correctly

### 05 Starter kit

The useful accessories kit provides support during daily wearing of compression stockings (e.g. gloves, lingerie bag, loyalty card)

### 06 Venotrain® Lotion\*

Developed in collaboration with Sebapharma; contains chamomile extract and is specially designed for people with sensitive skin; for use in the evening after taking off their stockings

\* only available in selected countries

# VenoTrain® glider

## COMPACT AID FOR PUTTING ON AND TAKING OFF STOCKINGS

The VenoTrain glider is a practical aid for putting on and taking off medical compression stockings. It can be prescribed by a physician.

The lightweight, hard-wearing material makes the stocking glide effortlessly over the arm or leg. It is suitable for compression stockings with an open or closed toe.

Saving space and time, it is the ideal companion when you're out and about.

### FEATURES

- Extremely lightweight, robust, and tear-resistant
- Easy to handle
- Ideal for travel



A video with donning instructions can be found at:  
[video.bauerfeind.com/venotrain-glider](https://video.bauerfeind.com/venotrain-glider)



PUTTING ON THE STOCKINGS: STEP 01



PUTTING ON THE STOCKINGS: STEP 02



PUTTING ON THE STOCKINGS: STEP 03



PUTTING ON THE STOCKINGS: STEP 04



PUTTING ON THE STOCKINGS: STEP 05



TAKING OFF THE STOCKINGS

# VenoTrain® glider plus

## ROBUST DORNING AID

The VenoTrain glider plus is a robust donning aid with an ergonomic design for all types of medical compression stockings. It can be prescribed by the physician.

Folding handles make it easy to put on the stockings, avoiding the need to bend down. It has been tested by an independent institute for ease of use, even for customers with restricted mobility.

### FEATURES

- Compact, ergonomic design
- Light metal frame with folding handles
- Easy to use without the need for any great effort



A video with donning instructions can be found at:  
[video.bauerfeind.com/venotrain-glider-plus](https://video.bauerfeind.com/venotrain-glider-plus)



PUTTING ON THE STOCKINGS: STEP 01



PUTTING ON THE STOCKINGS: STEP 02



PUTTING ON THE STOCKINGS: STEP 03



PUTTING ON THE STOCKINGS: STEP 04



PUTTING ON THE STOCKINGS: STEP 05



PUTTING ON THE STOCKINGS: STEP 06




# STANDARD PRODUCT RANGE

## ROUND KNIT PRODUCTS


DESIGN	SIZE	CCL			EDGING / SILICONE TOP BAND							TOE		
		1	2	3	Top band	Lace	Lace "sensitive"	Micro silicone dots	Silicone dots	Silicone dot pattern	Garter-belt	Closed	Open	
<b>VenoTrain</b> discrétion	AD	I-IV	●*	●**		●							●	●
	AG	I-IV	●*	●**			●		●				●	●
	AT	I-IV	●*	●**		●							●	●
<b>look</b>	AD	S-XL	●	●		●							●	
	AG	S-XL	●	●			●	●					●	
	AT	S-XL	●	●		●							●	
<b>micro</b>	AD	S-XL	●	●		●			●				●	●
	AG	S-XL	●	●		●		●	●	●			●	●
	AT	S-XL	●	●		●							●	●
	ATU	S-XL	●	●		●							●	●
<b>micro</b> Design Edition (Jive, Tango)	AD	S-XL	●	●		●							●	
	AG	S-XL	●	●			●			●			●	
	AT	S-XL	●	●		●							●	
<b>business</b>	AD	S-XL	●	●		●							●	
<b>soft</b>	AD	S-XL	●	●		●			●				●	●
	AG	S-XL	●	●		●			●	●			●	●
	AGHB	S-XL	●	●		●					●		●	●
	AT	S-XL	●	●		●							●	●
<b>pure</b>	AD	S-XL	●	●		●			●				●	●
	AG	S-XL	●	●		●			●				●	●
<b>soft S</b>	AD	XS-XL	●	●		●			●				●	●
	AG	XS-XL	●	●		●			●	●			●	●
<b>cocoon</b>	AD	S-XL		●		●							●	
<b>impuls</b>	AD	I-VI		●	●	●			●				●	●
	AG	I-VI		●	●	●			●				●	●
	AGHB	I-VI		●	●	●					●		●	●
<b>ulcertec</b> 39 moderate	AD	XS-XL	34-46 mmHg			●							Liner	Overstocking
<b>ulcertec</b> 46 strong	AD	XS-XL	34-46 mmHg			●							Liner	Overstocking

\* 10-15 mmHg    \*\* 15-20 mmHg


### DESIGN




**AD**  
Knee-high stocking




**AF**  
Mid-thigh stocking




**AG**  
Thigh-high stocking




**AGHB**  
Thigh-high stocking with garter belt




**AT**  
Pantyhose




**ATU**  
Maternity pantyhose




**ATH**  
Gent's pantyhose




**ATX**  
Control-top pantyhose



**ATE**  
Single-leg pantyhose



**Toe closed**



**Toe open**

# CUSTOM PRODUCT RANGE

## ROUND KNIT PRODUCTS

DESIGN		CCL			EDGING / SILICONE TOP BAND							TOE			
		1	2	3	Top band	Lace	Lace "sensitive"	Micro silicone dots	Silicone dots	Silicone dot pattern	Garter belt	Closed	Open	No foot	Stump
<b>VenoTrain</b>		1	2	3	Top band	Lace	Lace "sensitive"	Micro silicone dots	Silicone dots	Silicone dot pattern	Garter belt	Closed	Open	No foot	Stump
	<b>discrétion</b>	AD	•*	•**	•				•			•	•		
		AG	•*	•**			•			•			•	•	
	AT	•*	•**		•							•	•		
<b>look</b>	AD	•	•		•							•			
	AG	•	•		•		•	•		•		•			
	AT	•	•		•							•			
<b>micro</b>	AD	•	•		•				•			•	•		
	AF	•	•		•		•		•	•		•	•		
	AG	•	•		•		•		•	•		•	•		
	AT	•	•		•							•	•		
	ATU	•	•		•							•	•		
<b>micro Design Edition (Jive, Tango)</b>	AD	•	•		•							•			
	AG	•	•				•			•		•			
	AT	•	•		•							•			
<b>business</b>	AD	•	•		•				•			•			
<b>soft</b>	AD	•	•		•				•			•	•	•	•
	AF	•	•		•	•			•	•		•	•	•	•
	AG	•	•		•	•			•	•		•	•	•	•
	AGHB	•	•								•	•	•	•	•
	AT	•	•		•							•	•	•	•
	ATU	•	•		•							•	•	•	•
	ATH	•	•		•							•	•	•	•
	ATX	•	•		•							•	•	•	•
ATE	•	•		•							•	•	•	•	
<b>soft S</b>	AD	•	•		•				•			•	•		
	AG	•	•						•	•		•	•		
<b>impuls</b>	AD		•	•	•				•			•	•		
	AG		•	•	•				•			•	•		
	AGHB		•	•							•	•	•		
	AT		•	•	•							•	•		
	ATH		•	•	•							•	•		
	ATX		•	•	•							•	•		
<b>ulcertec 39 moderate</b>	AD	34–46 mmHg			•							Liner	Overstocking		
<b>ulcertec 46 strong</b>	AD	34–46 mmHg			•							Liner	Overstocking		

\* 10–15 mmHg \*\* 15–20 mmHg

Please refer to the current price list for information on standard and custom-made designs. The **custom-made lymphology product range** can be found in the current Lymphology Product Catalog (item number 271 4210 1001 001) as well as in the current price list. If you have any questions about our products, please contact our Customer Service.



# ADDITIONAL SERVICES ALTERATIONS | REPAIRS

## INDIVIDUAL CUSTOMER SERVICE

### ADDITIONAL SERVICES FOR NEW ORDERS

○ CUSTOM-MADE SOLUTIONS

	look	micro	micro DESIGN EDITION	business	soft	soft S	impuls
<b>Knitted toe area</b> Lengthen or shorten		○			○	○	○
<b>Toe</b> Place seam on the outside	○	○	○	○	○	○	○
<b>Toe</b> Wider, narrower, anatomical		○	○	○	○	○	○
<b>Panty area</b> Front/back – two seams	○	○			○*		○
<b>ATH</b> with longitudinal slit					○		○
<b>ATH</b> with transverse opening					○		○
<b>AT</b> with extra-wide gusset for the fuller figure	○	○			○		○
<b>Sew in a stoma opening</b> (only for repeat order)					○		

\* Please consult with Bauerfeind Customer Service.



## ALTERATIONS AND REPAIRS

● STANDARD

○ CUSTOM-MADE SOLUTIONS

		look	micro	micro DESIGN EDITION	business	soft	pure	soft S	cocoon	impuls	ulcertec
<b>Top of stocking</b>	Replace silicone top band	●○	●○	●○	○	●○	●○	●○		●○	
	Change top band to silicone top band	○	●○		○	●○		●○		●○	
	Make silicone top band narrower or wider	●○	●○	●○	○	●○	●○	●○		●○	
	Remove hip section and sew on silicone top band					●○				●○	
<b>AT</b>	Shorten panty area	●○	●○	●○		●○				○	
	Lengthen panty area					○*				○	
	Make gusset narrower (excluding ATX)	○	○*			○*				○	
	Make gusset wider (excluding ATX)	○	○*			○*				○	
	Extra-wide gusset for the fuller figure (excluding ATX)	○	○			○				○	
	Alter AT to AG with silicone top band	○	○	○		○				○	
<b>AGHB</b>	Replace hip strap					●○				●○	
	Shorten hip section					●○				●○	
	Shorten/Lengthen hip strap					○				○	
<b>Add a stoma opening or open the crotch afterwards (excluding ATX)</b>						○					
<b>Shorten stocking</b>	AD by up to 3 cm	●○	●○	●○		●○	●○	●○		●○	
	AF by up to 4 cm		○			○					
	AG by up to 5 cm	●○	●○	●○		●○	●○	●○		●○	
<b>Redye from light to dark**</b>						○		○		○	
<b>Toe alteration</b>	Alter closed toe to open toe		○*			○*	○*	○*		○*	
	Shorten open toe (from 3 cm)		○			○	○	○		○	
<b>Repairs, if possible</b>	Hole*	●○	●○		●○	●○		●○	●	●○	●○
	Run*	●○	●○		●○	●○		●○	●	●○	●○
	Panty area seam	●○	●○			●○				○	
	Loose thread*	●○	●○		●○	●○		●○	●	●○	●○

\* Please consult with Bauerfeind Customer Service.

\*\* Please do not try on stockings – re dyeing is no longer possible after contact with the skin.

Please refer to the current price list for information about designs, item numbers, and prices for standard and custom-made solutions.

The **custom-made lymphology product range** can be found in the current Lymphology Product Catalog (271 4210 1001 001) as well as in the current price list.

If you have any questions about our products, please contact Bauerfeind Customer Service.



211

Hospital


211

VenoTrain®	micro	soft	soft S	look	business	impuls
Ccl 1						
Ccl 2						
Ccl 3						
AD						
AF						
AG						
AGHB (right)						
AGHB (left)						
AT						
ATU						
ATH (without fly)						
ATX						
ATE						
Silicone dots (3 cm, only AD)						
Silicone dots						
Silicone dot pattern						
Micro silicone dots (only AG)						
Lace						
Lace sensitive						

Color:

(Mandatory information!)

open					
closed					



Account No.:

--	--	--	--	--	--

Order No.:

Name:

First Name:

Last Order No.:

<input type="checkbox"/>	Piece	<input type="checkbox"/>	Pair
--------------------------	-------	--------------------------	------

Panty Area (cG – cT) in cm:

front:                      back:

All sizes in cm. The purchaser is responsible for ensuring that all size information is correct. Custom-made solutions cannot be returned in exchange for a credit note or another item.

Special Requests:

Date

Stamp/Signature



# ORDER PADS

We have developed a range of order pads for you to order standard and custom-made solutions. These pads help you to order the required version of a particular VenoTrain product quickly and easily.

We are happy to send you these order pads free of charge.

They can also be downloaded at:

[www.bauerfeind.com/downloads](http://www.bauerfeind.com/downloads)

VenoTrain	Item number
Standard order pad	271 2020 0000 005
Order pad for custom-made solutions	271 2020 0001 003

## TERMS OF SALE AND DELIVERY OF BAUERFEIND AG (GENERAL BUSINESS TERMS) INTERNATIONAL

### Prices

are quoted net ex works in Euro currency, plus freight and taxes, if any.

### Dispatch

at buyer's cost and risk even if the buyer does not bear the transportation cost. Disruptions of operations affecting us or any of our suppliers, lack of suitable personnel, strike, lockout and force majeure are causes for extending delivery and relieve us of any delivery obligation for the time of the obstruction.

### Returns

outside the warranty terms are only accepted after prior advise and approval of Bauerfeind to the return, provided the delivery documents are submitted; it is supposed that the returned goods are in proper state. We reserve the right of rejecting the acceptance of any returned consignment. We charge an administrative fee of € 10 for every consignment returned outside the warranty period. The administrative fee for goods returned 1 month or more after the date of delivery is additionally 20 percent of the net order value, subject to a maximum of € 25 per article and piece; the fee is additionally 40 percent of the net order value of goods returned between 6 months and 12 months after the delivery date, subject to a maximum fee of € 50 for each article and piece. Goods returned outside the warranty, especially custom made products, and 12 months or more after the date of delivery will not be accepted back.

### Payment

to be made immediately after presentation of the invoice. First deliveries are only made against advance payment. If the customer is in default, any further delivery will only be made against advance payment. We charge interest on arrears at the rate of 8 percent above the applicable base rate from the time at which payment is delayed. Payments by check, draft, etc. are deemed to be made when the amount has finally been credited to our account.

### Retention of title

We retain the title to all goods delivered until the buyer has paid all amounts outstanding under our business relations with him. The retention of title also applies proportionately to goods which have been processed or otherwise modified. If the goods are sold all claims on third parties are deemed to have been assigned to us without proviso up to the total amount owing to us. Notwithstanding the assignment, the buyer is entitled to the collection of these amounts. Our right to collect the amounts is not affected by the preceding provision; however, we undertake to refrain from collection unless the buyer defaults on payment. If the buyer defaults on payment we can demand that the buyer disclose the assigned claims and the names of the debtors, provides all information required for collection of the debt, delivers all required documents and informs the debtors (third party) of the assignment. The buyer is not permitted to sell, hypothecate or assign the goods as security. Hypothecation, confiscation or any other action of constraint taken by a third party shall be communicated to us by the buyer

without delay and he shall provide all information and documents to enable us to defend our rights. We undertake to release, on the buyer's demand, any security which exceeds a secured claim that has not been satisfied by more than 20 percent.

### Complaints

Warranty claims due to obvious defects expire unless notified to us in writing without delay within the meaning of Article 377 HGB (German Commercial Code). For the rest, the time bar for warranty claims is 1 year after delivery of the goods purchased. If the goods are defective, we make replacement delivery or repair the goods at our choice. For reasons of hygiene, repairs can only be made to washed goods and goods which are of good hygienic condition. We will reserve the right of charging an adequate administration fee if unwashed goods are returned. If and to the extent to which we fulfil our duties of supplementary performance, in particular, the repair of defects or delivery of goods without defect, the customer cannot demand a price reduction or rescind the contract unless our attempts at supplementary performance have failed. Defects which are due to normal wear and tear, overstrain, misuse, neglect of care or non-compliance with the instructions for use are not covered by warranty.

### Liability

for damage is excluded unless we or our servants are found to have acted with gross negligence or intention. This limitation of liability does not apply to damage due to injury of life, body or health or the violation of a major contractual duty.

### Legal venue

The place of performance is the business seat of this company. The legal venue for all controversies from the business relations is Gera.

### Partial invalidity

If any provision in these Terms of Sale and Delivery should be or become invalid or unenforceable, the other provisions will not be affected by this. The ineffective or unenforceable provision will be substituted by a valid and enforceable provision which reflects the original intention of the parties as closely as possible.

### Data protection clause / call recording

For quality purposes, we randomly record calls between the buyer and our employees. The phone records exclusively serve the purpose of improving our performance and are deleted after evaluation. The data will not be disclosed to third parties. The buyer hereby consents to the records of telephone calls to that sole purpose. Consent is voluntary and may be revoked at any time in the future by the buyer.



# BETTER QUALITY OF LIFE FOR YOUR CUSTOMERS

Have you got any questions, suggestions, or requests?  
Please contact us:

## **BAUERFEIND AG**

Triebeser Straße 16  
07937 Zeulenroda-Triebes  
Germany  
**P** +49 (0) 36628 66-2000  
**F** +49 (0) 36628 66-2999  
**E** info@bauerfeind.com

## **AUSTRIA**

Bauerfeind Ges.m.b.H.  
Hainburger Straße 33  
1030 Wien  
**P** +43 (0) 800 4430-130  
**F** +43 (0) 800 4430-131  
**E** info@bauerfeind.at

## **BENELUX**

Bauerfeind Benelux B.V.  
Waarderveldweg 1  
2031 BK Haarlem

## **THE NETHERLANDS**

**P** +31 (0) 23 531-9427  
**F** +31 (0) 23 532-1970  
**E** info@bauerfeind.nl

## **BELGIUM**

**P** +32 (0) 2 527-4060  
**F** +32 (0) 2 792-5345  
**E** info@bauerfeind.be

## **BOSNIA AND HERZEGOVINA**

Bauerfeind d.o.o.  
Meše Selmovića 19  
71000 Sarajevo  
**P** +387 (0) 33 710-100  
**F** +387 (0) 33 619-422  
**E** info@bauerfeind.ba

## **CROATIA**

Bauerfeind d.o.o.  
Goleška 20  
10020 Zagreb  
**P** +385 (0) 1 6542-855  
**F** +385 (0) 1 6542-860  
**E** info@bauerfeind.hr

## **FRANCE**

Bauerfeind France S.A.R.L.  
B.P. 50258  
95957 Roissy CDG Cedex  
**P** +33 (0) 1 4863-2896  
**F** +33 (0) 1 4863-2963  
**E** info@bauerfeind.fr

## **ITALY**

Bauerfeind AG  
Filiale Italiana  
Piazza Don Enrico Mapelli 75  
20099 Sesto San Giovanni (MI)  
**P** +39 02 8977 6310  
**F** +39 02 8977 5900  
**E** info@bauerfeind.it

## **NORDIC**

Bauerfeind Nordic AB  
Kungsholmstorg 8  
11221 Stockholm  
**P** +46 (0) 774 100-020  
**F** +46 (0) 774 100-021  
**E** info@bauerfeind.se

## **REPUBLIC OF MACEDONIA**

Bauerfeind Dooel Skopje  
50 Divizija 24 a  
1000 Skopje  
**P** +389 (0) 2 3179-002  
**F** +389 (0) 2 3179-004  
**E** info@bauerfeind.mk

## **SERBIA**

Bauerfeind d.o.o.  
102 Omladinskih brigada  
11070 Novi Beograd  
**P** +381 (0) 11 2287-050  
**F** +381 (0) 11 2287-052  
**E** info@bauerfeind.rs

## **SINGAPORE**

Bauerfeind Singapore Pte Ltd.  
Blk 41 Cambridge Road #01-21  
Singapore 210041  
**P** +65 6396-3497  
**F** +65 6295-5062  
**E** info@bauerfeind.com.sg

## **SLOVENIA**

Bauerfeind d.o.o.  
Dolenjska cesta 242 b  
1000 Ljubljana  
**P** +386 (0) 1 4272-941  
**F** +386 (0) 1 4272-951  
**E** info@bauerfeind.si

## **SPAIN**

Bauerfeind Ibérica, S.A.  
Historiador Diago, 13  
46007 Valencia  
**P** +34 96 385-6633  
**F** +34 96 385-6699  
**E** info@bauerfeind.es

## **SWITZERLAND**

Bauerfeind AG  
Vorderi Böde 5  
5452 Oberrohrdorf  
**P** +41 (0) 56 485-8242  
**F** +41 (0) 56 485-8259  
**E** info@bauerfeind.ch

## **UNITED ARAB EMIRATES**

Bauerfeind Middle East FZ LLC  
Dubai Healthcare City  
Building 40, Office 510  
Dubai  
**P** +971 4 4335-684  
**F** +971 4 4370-344  
**E** info@bauerfeind.ae

## **UNITED KINGDOM**

Bauerfeind UK  
85 Tottenham Court Road  
London  
W1T 4TQ  
**P** +44 (0) 121 446-5353  
**F** +44 (0) 121 446-5454  
**E** info@bauerfeind.co.uk

## **USA**

Bauerfeind USA, Inc.  
1230 Peachtree St NE  
Suite 1000  
Atlanta, GA 30309  
**P** +1 800 423-3405  
**P** +1 404 201-7800  
**F** +1 404 201-7839  
**E** info@bauerfeindusa.com